



Bore fitting the Stryker handprobe

BPEI's
Stryker
handprobe



1st generation MOOKP Micro-Mill
demonstrated at the 8th KPro meeting,
May 12, 2012



2nd generation MOOKP Micro-Mill

2nd generation MOOKP micro-Mill with Stryker CORE Impaction Drill



Sterilization: EtO gas and hydrogen peroxide gas plasma (e.g., Sterrad 100 S)

Probe: Stryker CORE Impaction Drill. For other drill brands, the micro-Mill needs modifications.

For more information, contact abernal@bioniko.com